


An exploded view of a mechanical device, possibly a medical instrument, shown in a light grey color. The device is composed of numerous parts, including a large L-shaped frame, a circular component with a lens or sensor, and various smaller components like screws, bolts, and tubes. The background is a vibrant green and blue gradient with abstract circular patterns and a grid of dots.

BIOVOLUME®

A GLOBALLY CONNECTED AI
IMAGING AND ANALYSIS PLATFORM
REVOLUTIONISING THE COLLECTION
OF SUBCUTANEOUS TUMOUR DATA
IN PRE-CLINICAL ONCOLOGY.

MORE THAN 50% OF PRE-CLINICAL STUDIES CANNOT BE REPRODUCED



COSTING THE
PRE-CLINICAL
INDUSTRY
\$28.8BN
ANNUALLY.

PLOS Biology, The Economics
of Reproducibility in
Preclinical Research

**EFFECTS OF ANTI-CANCER
DRUGS CAN BE EASILY
MISSED**

[WITH CALIPERS] AS TUMORS WILL TEND
TOWARDS BEING DETERMINED WITH A GREATER
BIAS AS THEY GROW LARGER.

Ishimori et al. Acad Radiol 2005; 12:776-781

CREATING A COMMERCIAL CHALLENGE




5 IN 5,000

DRUGS THAT BEGIN PRE-CLINICAL TESTING MAKE IT TO HUMAN TESTING

1 IN 5

WILL NEVER BE APPROVED FOR HUMAN USE



Hollingshead. JNCI. Vol. 100, Issue 21, November 5, 2008

Scientific American – Cost to Develop New Pharmaceutical Drug Now Exceeds \$2.5B. – Rick Mullin.


BIOVOLUME PROCESS

SCANNER,
LAPTOP



AZURE CLOUD PLATFORM

3D STEREO



ANIMAL ID	SCAN DATE/TIME	GROUP NAME	LENGTH (mm)	WIDTH (mm)	HEIGHT (mm)	VOLUME (mm ³)	OPERATOR (user)	MANUAL SEGMENTATION	PROCESSING	EXCLUDE
5001	22/02/2019	AU484	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
5002	22/02/2019	AU485	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input checked="" type="checkbox"/>
5003	13/02/2019	AU481	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
112	02/02/2019	AU482	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
5004	22/01/2019	AU483	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
34	13/01/2019	AU485	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
340	22/02/2019	AU486	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
341	22/02/2019	AU487	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
431	22/02/2019	AU488	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>
50	22/02/2019	AU489	7.04	5.1	101	6006.125	Andy Smith	Yes	✓	<input type="checkbox"/>

THERMOGRAPHY

3D & THERMAL
CAPTURE

RECONSTRUCT

SEGMENT

ANALYSE

IMPORT
EXPORT
API

STATISTICS &
MACHINE
LEARNING
ON AZURE


GREATER TRACEABILITY - THROUGH A CLOUD PLATFORM

- + THERMAL CAPTURE
- + MULTI-TENANTED CLOUD PROCESSING
- + BIG DATA SET
- + BROWSER ENABLED
- + CONSOLIDATED STUDY DATA & ANALYTICS DASHBOARD



Additional Thermal Capture to aid accuracy and precision

MUCH MORE THAN CALLIPERS



GLOBAL DATASET NOW IN EXCESS OF 10,000 SCANS

31ST May

Total scans: 6,532

8TH August

Total scans: 8,561

13th September



Total scans:

10,806


Scans #	Tenants #	Strains #	Users #
10,806	9	12	50

Tenant ID All	Bin size 100
------------------	-----------------

Scans per volume
Volume (bin)



Animal Strain Name			Scans / tumour / day		
Null	Fisher Rat	SCID	1	2	3
ATH	FVB/N		4	5	6
Balb/C	NSG		7	8	9
BL6	NSG-A2		10	11	12
C57BL6	Nude		13	14	15
C57Bl/6	R2G2		16	17	18
			19	20	21
			22	23	24



IMPROVING THE SCIENCE THROUGH DIGITISATION AND DATA

BASED ON CURRENT FORECASTS WE WILL HAVE




IN OPERATION MID-2020

NUMBER OF UNITS	200	
EACH UNIT TAKES	80	A DAY
EACH UNIT USED	5	DAYS A WEEK
UNIT USED OVER	12	MONTHS

TOTAL DATA COLLECTED

960,000

TRANSFORMING THE SCIENCE





A fully traceable cloud connect platform is designed to transform the way drug efficacy is monitored and understood

CONTINUOUS LEARNING

ARTIFICIAL INTELLIGENCE EFFORT

USES OF AI AT FUEL3D



PLOS ONE

RESEARCH ARTICLE

An innovative non-invasive technique for subcutaneous tumour measurements

Juan Delgado-Santamaría^{1,2,3*}, Gaelle Ehrhardt^{4,5}, Martin Paszkowski¹, Sam Hackett¹, Andrew Smith¹, Wojciech Winiarski⁶, James Klotzow⁷, Adalberto Zúñiga⁸, Anna Chłobudzka⁹, Leonardo Rubio-Naveira¹⁰, Anwar Faheem¹¹, Gema Wilson¹²

1 Fuel3D, Oxford Science Park, Oxford, England, United Kingdom, **2** AstraZeneca MSD Biotech Unit, Bioregry Sciences, Cambridge Science Park, Cambridge, England, United Kingdom, **3** AstraZeneca MSD Biotech Unit, Oncology, Addenbrooke Park, Macclesfield, England, United Kingdom, **4** Current address: Analytic Department, Wroclawskie Centrum Badań, Wrocław, Poland, **5** Current address: Institute for Mathematical Innovation, University of Bath, Bath, England, United Kingdom, **6** Current address: DENSO, Lindeu (Bismarck), Germany, **7** Current address: CEM, AstraZeneca, England, United Kingdom, **8** Current address: Shearwater, London, England, United Kingdom, **9** Current address: AstraZeneca, England, United Kingdom, **10** Current address: AstraZeneca, England, United Kingdom, **11** Current address: AstraZeneca, England, United Kingdom, **12** Current address: AstraZeneca, England, United Kingdom

* jdelgado@fuel3d.com

OPEN ACCESS

Citation: Delgado-Santamaría J, Ehrhardt G, Paszkowski M, Hackett S, Smith A, Winiarski W, et al. (2019) An innovative non-invasive technique for subcutaneous tumour measurements. *PLoS ONE* 14(10): e0216690. <https://doi.org/10.1371/journal.pone.0216690>

Editor: Matthew Ryan, Cardiff University, UNITED KINGDOM

Received: April 14, 2019

Accepted: September 10, 2019

Published: October 14, 2019

Copyright: © 2019 Delgado-Santamaría et al. This is an open access article distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Data Availability Statement: All student data is available on our website www.fuel3d.com for the real world data. All data files are available in the study setup, data collection and analysis, desktop or tablet, or preparation of

Abstract

In oncological drug development, animal studies continue to play a central role in which the volume of subcutaneous tumours is monitored to assess the efficacy of new drugs. The tumour volume is estimated by taking the volume to be that of a regular spheroid with the same dimensions. However, this method is subjective, insufficiently traceable, and is subject to error in the accuracy of volume estimates as tumours are frequently irregular.

Methods & results

This paper reviews the standard technique for tumour volume assessment, calliper measurements, by considering a statistical review of a large dataset consisting of 2,500 tumour volume measurements from 1,600 mice by multiple operators across 6 mouse strains and 20 tumour models. Additionally, we explore the impact of six different tumour morphologies on volume estimation and the detection of treatment effects using a computational tumour growth model. Finally, we propose an alternative method to callipers for estimating volume—BioVolume™, a 3D scanning technique. BioVolume simultaneously captures both stereo (RGB) (Red, Green and Blue) images from different light sources and infrared thermal images of the tumour in under a second. It then detects the tumour region automatically and estimates the tumour volume in under a minute. Furthermore, images can be processed in parallel within the cloud and as the time required to process multiple images is similar to that required for a single image. We present data of a pre-production unit test consisting of 297 scans from over 120 mice collected by four different operators.

Conclusion

This work demonstrates that it is possible to record tumour measurements in a rapid, minimally invasive, morphology-independent way, and with less human bias compared to

PLOS ONE | <https://doi.org/10.1371/journal.pone.0216690> | October 14, 2019 | 1/14

PLOS ONE

Subcutaneous tumour scanner




Fig 1. Repetitive snap capture set up of the BioVolume unit including computer monitor and desktop (left), along image of a white XCD mouse being processed by the scanner (right).

<https://doi.org/10.1371/journal.pone.0216690.g001>

Abstract

region is exposed to the device operator (see Fig 1, left). Then, acquisition is triggered on either a cable-connected laptop or by using a button on the device itself (see Fig 1). There is no requirement to anaesthetise the animal. Acquisition takes around 0.2s, and rendering occurs in the cloud. Rendering a full image, segmentation, and measurement extraction requires approximately 2s, depending on internet speed and these can be parallelised. The acquired measurements are then displayed to the operator immediately.

The BioVolume unit consists of a stereo system with two RGB cameras, three white light flashes, and an infrared thermal camera. Upon activation, the unit collects 6 photographic (RGB) images and a single thermal frame (see Fig 2, left), for a primary flash and thermal image example). The flash of a stereo software is linked to the work with a beta version named v0.1_0415. The RGB images are reconstructed in a surface by means of a binocular stereo-process, acquiring both depth and RGB maps (Fig 2, right). The thermal frame is mapped onto the 3D reconstruction using a conventional affine transformation based on a prior positional calibration of the RGB and thermal camera. The segmentation of the tumour happens on the thermal map, which is then projected onto the depth map. The height is then obtained by fitting a plane to the back of the mouse using an optimisation algorithm. More details are provided in [S1 Appendix](https://doi.org/10.1371/journal.pone.0216690).

Volume calculation. We compare two formulas (1,2) for the estimation of tumour volume:

- Spheroid formula (BioVolume & calliper):
$$V_{\text{sph}} = \frac{4}{3} \pi \cdot \text{length} \cdot \text{width}^2 \quad (1)$$
- Cylindrical volume (BioVolume):
$$V_{\text{cyl}} = \text{Area} \cdot \text{height} \quad (2)$$

Computational to mouse growth model. The cellular automation model consisted of a rule-based model operating on two simulated cell populations growing on a 3D lattice. The

PLOS ONE | <https://doi.org/10.1371/journal.pone.0216690> | October 14, 2019 | 2/14

PLOS ONE

Subcutaneous tumour scanner



Fig 2. Summary of tumour acquisition and XCD in volume from calliper. Panel A: Quantification of tumour volume from a large set of XCD values in mouse of different sizes. Panel B: Values plotted on the percentage number of observations. Panel C: Values plotted on the percentage number of observations. Panel D: Values plotted on the percentage number of observations.

<https://doi.org/10.1371/journal.pone.0216690.g002>

at least half the reported volume. In contrast, only 38.8% of weights exhibited relative errors of less than 50% and only 17.4% of comparisons returned errors of less than 20%.

Cellular automation model

Tumour morphology. To assess the impact of using calliper measurements to estimate tumour volume we developed a Cellular Automata (CA) model of tumour growth and its treatment (see Fig 3).

We simulated the growth of both control and treated tumours for each morphology described in Fig 4. Treatment commenced on day 15 and lasted for 10 days. The tumours were subject to standard calliper (CC) measurements which were then used to estimate the spheroidal volume of the tumour. This volume was then compared to the actual tumour volume (computed as

PLOS ONE | <https://doi.org/10.1371/journal.pone.0216690> | October 14, 2019 | 3/14

<https://journals.plos.org/plosone/article/file?id=10.1371/journal.pone.0216690&type=printable>

MORE DATA, BETTER SCIENCE: A GREATER ADVANTAGE



REPRODUCIBLE RESULTS



GREATER SCIENTIFIC CONFIDENCE



BETTER TRACEABILITY

INTRODUCING NEW BIOMARKERS:



TUMOUR CONDITION



VASCULAR PHYSIOLOGICAL MARKERS



An exploded view of a mechanical assembly, possibly a medical device component, rendered in a light beige color. The assembly is shown in a disassembled state, with various parts like screws, bolts, washers, and structural plates scattered around the main components. The background is a gradient of green and blue with abstract circular patterns.

BIOVOLUME[®]

A GLOBALLY CONNECTED AI
IMAGING AND ANALYSIS PLATFORM
REVOLUTIONISING THE COLLECTION
OF SUBCUTANEOUS TUMOUR DATA
IN PRE-CLINICAL ONCOLOGY.