

# Capsuled

## True data-driven patient support

A personally customized medication adherence solution supporting a specific patient journey for any medication modality.

Designed to create an accessible gateway to all digital information for the patient, provided at the right time, along with timed encouragement messaging.

Capsuled makes use of auditory and visual alerts to remind a patient to take medication at the designated time, integrates educational videos and closes the patient-caregiver-Patient Support Program loop to enhance care outcomes.



#### How does Capsuled work?



Upon medication extraction, the unique med ID will be sent to Vaica's cloud

#### 01 💂

Patient Support Program remotely sets specific device settings via Vaica's portal

patient settings

02



If a patient doesn't take
the med on time,
an e-mail / SMS will be sent
automatically to the patient,
PSP's and/or caregivers for
further intervention



03 😉

Slightly before scheduled intake time, the screen shows a "Take med" color escalating notification

Once activated, Capsuled

automatically downloads

**i** 0

If the patient takes the med on time, it will be logged on the adherence report



When the scheduled time arrives an audible alert is sounded

### Capsuled Features:

• Real time reminders enabling patients, PSP and caregivers real time interventions

• 4G connectivity cloud services based on an internal cellular modem

• All medication modalities inhalers, injectable, orals

• Interactive questionnaires PSP feedback, patient engagement

• Timed encouragement messaging full support system

Integrated educational videos \*optional personalized technology

Locking Mechanism \*optional

• Built-in hub \*optional for vital sign monitoring devices

• Built-in NFC reader \*optional unique med ID identification

**Dimensions:** Customizable | **Temperature:** -10°C to 45°C (-50°F to 110°F)

Main Power: 120-240 VAC, 50/60 Hz | Backup Power: Li-ion Battery, Up to 10 day backup

Communication: 3G / 4G, BLE, NFC

